

cobas e 411 analyzer

Peace of mind in immunology





Sample loading in STAT mode



Ready to load epack

Successor to the Elecsys[®] 2010, one of the most successful immunoanalyzers ever developed, the **cobas e** 411 analyzer is a perfect complement to the **cobas c** 311 analyzer.

- The cobas e reagent packs offer best-in-class handling convenience and reagent efficiences for immunochemistry testing with liquid ready-to-use reagents and 60 day average onboard reagent stability
- 9 minute STAT applications for Troponin, CK-MB Myoglobin, hCG and PTH provide unmatched response times
- As the heart of the cobas e 411 analyzer, ECL technology delivers the unique combination of high analytical sensitivity, wide measuring ranges and 9 minute STAT applications for unparalleled performance in immunochemistry

_ 2 of 9





- CD-ROM, USB & printer ports enable efficient data management
- 30 continuous random access sample positions with STAT priority
- Also available rack version with 75 sample positions with STAT priority

- cobas* touchscreen user interface is standardized to other cobas systems
- 18 channels enable use of up to 18 parameters simultaneously
- 6 Easy removal of disposables via solid waste container

- Disposable tips and cups eliminate carryover
- Measuring cell provides unmatched detection using ECL technology
- Consumables and liquid waste are easily accessible

Full workstation consolidation

- The cobas e 411 analyzer is designed to deliver accurate results on account of proven ECL technology
- It consolidates routine endocrinology with cardiac, cancer, bone and infectious disease markers

Serum work area integration

- Combined with the cobas c 311 analyzer, the system processes up to 95% of your serum work area workload
- Work area management solutions integrate sample and data management across systems

Seamless laboratory networking

- As part of the cobas modular platform, the cobas e 411 analyzer offers consistent results with other cobas e systems
- Universal reagent carriers and analytical components enable flexible network solutions








Enhanced data review screen

Convenience

cobas e packs for convenient handling

- · cobas e reagent packs are liquid, ready-to-use and simplify handling
- · All-in-one cassette packaging for each parameter combined with calibrators streamlines logistics
- · Automatic pack opening/closing ensures longer onboard stability
- · Programming-by-loading concept ensures seamless, precise data handling
- · Data review screen provides quick result traceability






Elecsys assay performance in external ring-trials

Performance

9 min. applications for speed and quality

- · Over 96 assays offer comprehensive coverage more than 10 indication areas
- · High onboard stability and shelf-life allow continuous availability of both routine and esoteric parameters
- · Disposable tips and assay cups for reaction eliminate risk of carryover
- · Integrated safety features ensure result integrity





Confidence

ECL technology: The gold standard

· High analytical sensitivity enables large measuring ranges and low sample volume

Measuring ranges: Estradiol 5-3,000 pg/mL; Progesterone 0.03-60.0 ng/mL;

Total PSA 0.002-100 ng/mL β-hCG requires 10 μL;

Sample Volume: 3rd Gen TSH requires 50µL

- · Activation with voltage for a controlled reaction allows high precision to avoid unnecessary repeats and reruns
- · Short incubation time enables the timely availability of innovative applications



cobase411 Assay Menu⁺

analyzer

Thyroid function	Tumor markers	Infectious diseases
Anti-TG	AFP	Anti-HAV
Anti-TPO	CA 125 II	Anti-HAV IgM
Anti-TSHR	CA 15-3	Anti-HBc
Calcitonin	CA 19-9	Anti-HBc IgM*
FT3 III	CA 72-4	Anti-HBe
FT4 II	CEA	Anti-HBs II*
T3	CYFRA 21-1	Anti-HCV II
T4	free PSA	HBeAg
T-Uptake	NSE	HBsAg II
Tg II	total PSA	HBsAg confirmatory*
Tg Confirmatory Test	Pro GRP	HbsAg II quant
TSH	He4	CMV IgG Avidity
	SCC*	CMV IgM
		CMV IgG
Fertility / hormones	Cardiac	HSV-1 IgG
ACTH	CK-MB (mass)	HSV-2 IgG
C peptide	CK-MB (mass) STAT	HIV Antigen confirmatory*
Cortisol II	Digoxin	HIV Antigen
DHEA-S	Digitoxin	Syphilis
Estradiol III	Myoglobin	HIV combi PT
FSH	Myoglobin STAT	Rubella IgG
hCG+β	NT-proBNP	Rubella IgM
hCG STAT	Troponin T	Toxo IgG Avidity
Insulin	Troponin T hs	Toxo IgM
LH	Troponin T hs STAT	Toxo IgG
Progesterone		
Prolactin		
SHBG		
Testosterone		
AMH		
Maternal care	Anemia	Bone markers
free β hCG	Ferritin	β-CrossLaps
PAPP-A	Folate III	Vitamin D Total
PIGF	RBC Folate	Intact PTH
sFLT-1	Vitamin B12	N-MID Osteocalcin
		PINP
		PTH STAT
Rheumatoid Arthritis	Critical care	
Anti-CCP	IL6	
	Procalcitonin	Others
	S100	IgE
Immuno Suppressive Drug		hGH
Tacrolimus		
Cyclosporine		
Sirolimus		
Everolimus		

cobas e 411 analyzer

Technical Specification



System	Fully automated, continuous loading and random access system for immunoassay analysis based on ECL technology. Available as a disk and rack system. System components: Analytical module including touchscreen PC		
Test Throughput	Up to 86 tests/hr		
Reagent concept	cobas e reagent pack: Ready-to-use reagents Onboard reagent capacity: 18 channels/reagent slots for up to 18 different assays Evaporation protection: Automatic closure of reagent lids		
Programmable parameters	Max 60 assays definable via 2D-barcode including all necessary information (programming by loading)		
Sample types	Serum, Plasma, Urine and others		
Onboard sample capacity	Load/unload capacity: Rack types: STAT handling:	Disk: 30 samples Rack: 75 samples on 15 RD standard 5 position rack Routine, STAT, Control, Calibrator Disk: Any unoccupied position on the sample disk Rack: Dedicated STAT port on rack feeder	
Sample container types	Primary tubes: Sample cup: Cups on tube	5–10 mL; 16x100, 16x75, 13x100, 13x75 mm 2.5 mL	
Sample volume	10 to 50 μL per test, depending on assay protocol		
Sample probe	Disposable tips for individual sample to avoid carry over; Sample clot detection; Liquid level detection		
Sample barcode types	PDF417, Code 128, Cod	PDF417, Code 128, Codabar (NW 7), Interleaved 2 of 5, Code 39	
System interfaces	Touchscreen PC, RS 232 serial interface, bi-directional Standard PC ports (USB, ethernet, serial etc) for other communication devices		
Sample database	2,000 results for routine, STAT and controls		
Test methods	Pre-defined assay prot	Pre-defined assay protocols (sandwich, competitive, titration)	
Calibration methods	2-point calibration per lot or per cobas e reagent pack		
QC methods	Individual QC; Cumulative QC; Preventive QC after calibration		
Rerun/repeat/reflex	Rerun, repeat and reflex supported		
Electrical requirements	Power requirements: Power consumption:	100–120 VAC 50/60 Hz single phase or 200–240 VAC 50/60 Hz single phase 1,000 VA	
Water/waste requirements	Water supply: Water quality: Liquid waste: Water consumption:	3l container ≤10µS/cm or ≥0.1 mega Ohm, bacteria-free 4l container with overflow detection or direct drain Approx. 3l for 250 tests/approx. 12 mL/cycle	
Regulatory compliance	GS, CE, UL, C-UL, CB-	GS, CE, UL, C-UL, CB-report and certificate	
Operating conditions	Ambient temperature: Ambient humidity: Noise output:	8 to 32°C (64.4°F to 89.6°F) 20% to 80% 60 dbA (stand-by mode), 63 dbA (avg. during operation)	
Physical dimensions	Width: Depth: Height: Weight:	120 to 170 cm (disk/rack) 73 to 95 cm (disk/rack) 56 cm (w/o PC unit) Approx. 170 kg (disk) and 210 kg (rack)	

COBAS, COBAS C, COBAS E, LIFE NEEDS ANSWERS and ELECSYS are trademarks of Roche.

Other brand or product names are trademarks of their respective holders.
© 2016 Roche Diagnostics

Roche Diagnostics India Pvt. Ltd.